## Disassembly and Recycling instructions

Knú Comfort encourages you to recycle whenever possible



**Pull Out Sleepers:** including; EV Collection, Florin, Harmony, Rêma Series, **Tranquility\*** (\*discontinued product)

#### What to expect:

**Disassembly Time:** 15-20 minutes per chair (General)

(allow additional 20-30 minutes for detailed component disassembly)

#### **Tools Needed**

Personal Protective Equipment (Safety Glasses) Torx Bit (T27) Socket/Ratchet or Wrench (7/16) Screwdrivers (Phillips #2 - #3) Utility Knife

#### **Material Content:**

Wood Metal (3 steel) Poly Foam Textile (fabric/vinyl) Polyethylene Plastic

#### **Component Content:**

<u>Upholstered Arms</u>: (plywood frame w/tnuts and staples, poly foam, textile w/staples (fabric/vinyl), arm caps (solid wood or urethane/steel) <u>Upholstered Back</u> (plywood frame w/t-nuts, poly foam, textile w/staples (fabric/vinyl), mounting brackets (steel)

<u>Upholstered Seat/Seat Extension:</u> (plywood w/t-nuts (steel), poly foam, textile w/staples (fabric/vinyl)

Pull-Out Mechanism: ( steel, plastic)

Pull-out Drawer Assembly: (plywood frame w/staple, drawer glides ( steel)

Hardware: (bolts/screws ( steel)

# **Pull Out Sleeper:**

**Mechanism Removal: (typical)** 

| • | Use screwdriver | (phillips) to | remove | (4) mounting | bolts | and r | emove | the t | front |
|---|-----------------|---------------|--------|--------------|-------|-------|-------|-------|-------|
|   | pull out panel  |               |        |              |       |       |       |       |       |



To allow access to the mechanism, place unit on it's back



 Lift/Pull out the seat/seat extension and allow seats to tilt and rest against back cushion





# Disassembly and Recycling instructions

Knú Comfort encourages you to recycle whenever possible





# Pull Out Sleeper: (con't)

| 1 1                             | m out siceper.  |   |  |  |  |  |
|---------------------------------|---|---|--|--|--|--|
| •                               | Use screw driver (phillips) to remove screw holding each plastic stop in place and remove from mech                                     |   |  |  |  |  |
|                                 |   |   |  |  |  |  |
| •                               | Lift/pull the seat/seat extension and mechanism as one unit   |   |  |  |  |  |
|                                 |   | CANCEL STORES   |  |  |  |  |
| •                               | Use screwdriver (phillips) to remove (4) mounting bolts and remove the seat extension   |   |  |  |  |  |
|                                 |   |   |  |  |  |  |
| •                               | Use screwdriver (phillips) to remove (4) mounting bolts and remove the pull out mechanism and recycle (6 steel)                         |   |  |  |  |  |
|                                 |   |   |  |  |  |  |
| •                               | Use screwdriver (phillips) to remove (3) mounting screws from each side and recycle (3 steel)   |   |  |  |  |  |
| Outside Beek Demovale (typical) |   |   |  |  |  |  |
|                                 | Place the unit on end (arm), Use screwdriver (phillips) to remove (3) mounting bolts from each arm and remove the outside back assembly |   |  |  |  |  |
| _                               |   | CONTRACTOR OF THE PARTY OF THE |  |  |  |  |

## **Inside Back Removal: (typical)**

• Remove inside back by sliding the back frame off the pivot mounting plate.



## Disassembly and Recycling instructions

Knú Comfort encourages you to recycle whenever possible





# Pull Out Sleeper: (con't)

### **Caster Removal: (typical)**

• Use socket/wrench (7/16") to remove (4) mounting bolts from each caster



### Leg Removal: (typical)

• Use screwdriver (phillips) to remove (2) mounting bolts from each leg



