



Hydraulic Lift Exam Recliner/Rema Style » **R7H67RU**.



Hydraulic Pedal (Left Side Facing) & Ballistic Nylon Protection



Quad Locking Foot Pedal & 3" Dual Wheel Caster

Hydraulic Lift Exam Recliner

The Hydraulic Lift Exam Recliner allows the caregiver to intuitively raise and lower the recliner from seated height, to exam height, and stop anywhere in between, while allowing the patient to operate the recliner manually. It comes standard with many important features, and is offered with extensive option combinations to provide the end user a chance to specify the perfect product for their application. Shown here is the Hydraulic Foot Pedal which is standard on each side of the recliner. There is a ballistic nylon material sewn into the bottom of the arm panel to help protect the chair material from extensive usage of the lift, Trendelenburg, or Quad Locking Caster pedals.

The Quad Locking Caster System is an important feature on a Lift Exam Recliner to provide stability when raising a patient up for examination, and is a standard option on all of our Lift Exam Recliners. The Quad Locking Caster System features 3" dual wheel rubber casters, which provide ultimate traction and stability. The Quad Locking Foot Pedal shown in the above image is standard on both sides of the recliner, in the front for easy access to lock all wheels with one action.