



Power Lift Exam Recliner/Rema Style » **R7P18J**.



Standard Trendelenberg Pedal



Power Lift Exam Activation Pedal

Power Lift Exam Recliner

The Power Lift Exam Recliner offers the same patient lifting functionality of the Hydraulic model, with the addition of the Power Motor and Foot Pedal Activation in lieu of the Hydraulic Pedal for a smooth and seamless travel range. It also comes standard with all the feature and option combinations as our Hydraulic Lift Exam Recliner. The Power Lift Exam Recliner is well suited for Hybrid Exam/Triage/Isolation Room applications where you can get all the functionality of an exam table, while taking up less space and providing recliner functionality when an Exam Table Treatment position is not needed.

Power Exam Lift Activation Pedal

The Power Exam Lift Activation Pedal is marked with a "+" and a "-" which will raise and lower the recliner when activated by stepping on the pedal. The Activation Pedals are located on both sides of the recliner. The Power Lift Motor features an anti-crush emergency stop, which will halt the motor if an object comes in contact with the chair while it is being lowered back down. The Power Exam Recliner also features an auto-off cycle function, to protect the motor.