

#### **Series Overview**

The Durable® Power Series combines over a half century of recliner mastery with leading, proven technology to create something completely new to the healthcare space—fully adjustable infinite position power recliners featuring weightless zero gravity capability. The intuitive control wand is preset with easy access buttons for the most common positions, while individual movement controls allow users to fine tune the footrest, back, and seat angle to their specific level of comfort. These features are found across healthcare recliners, exam lift recliners, and healthcare recliners with lift assist functionality.

#### **General Product Information**

- 650 lb Static / 500 lb Dynamic Weight Rating
- Designed for 24/7 use healthcare environments
- Styles include recliners, recliners with lift assist function, & exam lift treatment recliners
- Power operated recline functionality via control wand
- Footrest, back, and seat pitch can be controlled independently and precisely
- Control wand preset buttons include: Trendelenburg, Sleep, Zero Gravity, and Home (base position)

- Quick removable back & lift up seat for ease of cleaning
- True flatbed sleep position
- Standard with Chaise Continuous Leg Support feature
- Sinuous spring seat and back
- Seat clean-out
- Mid-Ottoman Calf Support
- Infinite Position Back reclines to infinite angle positions
- All interior plywood California Carb Compliant/Phase II
- Field replaceable components

#### Options vary by model, see specific model page for list of available options

Refer to pg 99 for full option details and images

# To Order, Specify

- Model Number
- Textile
- Finish
- Options







# **ZG** Power Healthcare Recliner

# **Select Aesthetic Style** (Denotes back & arm design)








(Koncert style pictured)

**Functional Style** (Aesthetic style code precedes selected model number below, e.g., KCZP17)

#### CZP17 ZG Power Healthcare Recliner

СОМ	Dillon/Carson	$A_5$	B <sub>6</sub>	C <sub>7</sub>	D <sub>8</sub>	E9	F <sub>10</sub>	$G_{11}$	H <sub>12</sub>	I <sub>13</sub>
9199	9479	9739	9999	10259	10529	10919	11319	11839	12369	13029

#### **Standard Features**





ARM CAP OPTIONS



Recline Functions Operated by Intuitive Control Wand

See pg 2 for finishes



4" Central Locking Caster System



Weightless Zero Gravity Position via Wand Preset

PATIENT WELLNESS & COMFORT OPTIONS

Di	m	er	ารเ	0	n	S

	Overall —			_	– Seat –			Kick Out			
	W	D	Н	W	D	Н	Back H	Arm H	Seat D	Recline D	from Wall
RCZP17   KCZP17	33	37	46.5	24	21	22	27.5	26.5	40.5	68	15.5
KMCZP17	33	37	46.5	24	21	22	27.5	26.5	40.5	68	15.5



# **Available Options**

U	Urethane	169	ABP	Adjustable Back Pad	
J1	Solid Surface - Group 1	619		Textile Grades COM-B	459
J2	Solid Surface - Group 2	819		Textile Grades C-I	569
F	Fully Upholstered Arms	Standard	AHP	Adjustable Head Rest Pillow	339
			DP	Dual IV Poles & Brackets <sup>1</sup>	498
			FBH	Foley Bag Holder	89
	CASTER OPTIONS		FT	Foot Tray	279
7.02	3" Decal Miles of Control Londring Sections	Na Chausa	IV	IV Pole & Bracket <sup>1</sup>	249
Z-Q3	3" - Dual Wheel Central Locking System	No Charge	LC	LiquiCell	229
Z-Q4	<b>O</b> ,	Standard	LCKFR	Orthopedic Locking Footrest <sup>6</sup>	Standard
Z-Q5	5" - Central Locking System	119	MBB	Moisture Barrier (Back)	59
			MBS	Moisture Barrier (Seat)	59
			MR	Monitor Arm (right side facing) <sup>1,2</sup>	849
	LIEAT/MASSACE ORTIONS (IEC COCOL 1	CERTIFIED)	MTR	Monitor Arm with Tray (right side facing) <sup>1,2</sup>	949
	HEAT/MASSAGE OPTIONS (IEC 60601-1	CERTIFIED)	PR1	Paper Roll Holder <sup>3</sup>	229
	Wand Controls		TE	Trendelenburg	Standard
HW	Heat	1139	USB	USB Power Outlet (IEC 60601-1 Certified)	459
HMW	Heat & Massage	1369	WR	Wrist Restraint Mounts (pair)	229
	In-Arm Controls			TRAY TABLE OPTIONS	
Н	Heat	849		TRAT TABLE OF HONS	
HM	Heat & Massage	1119	TSR1	Adjustable (Right) Side Tray <sup>1,2,5</sup>	399
			TD3	Dual Adjustable Side Trays (TSR1 & TSL1) <sup>1,5</sup>	798
	ADDITIONAL FEATURES & OPTIONS		LFSTR	Large Folding (Right) Side Tray <sup>4</sup>	249
BBU	Battery Back Up	649	LFSTRCH	Large Folding Cup Holder (Right) Side Tray	319
PB	Push Bar	Standard	PLFSTR	Padded Large Folding (Right) Side Tray	379
EG	Protective Edge Guard	Standard	SFSTR	Small Folding (Right) Side Tray <sup>4</sup>	249
QR	Additional Quick Release Accessory Brack				
CAL	CAL TB133	See pg 9	*D.	efer to pg 99 for full option details and images*	
Z	Multi Cover Use Pricing for Hi	, 0	l Re	ejer to pg 99 jor juli option details and images."	

<sup>&</sup>lt;sup>1</sup>Includes QR bracket connection installed on recliner

<sup>&</sup>lt;sup>6</sup>Position is locked in place by motors and can only be adjusted by the control wand.







 $<sup>^2</sup>Standard$  on right side (facing). Replace "R" with "L" for left side (facing) option where available.

 $<sup>^3</sup>$ Takes place of Push Bar. Paper sold separately.

<sup>&</sup>lt;sup>4</sup>Contact customer service for solid surface pricing

 $<sup>^5\</sup>mbox{Over the lap tray function recommended in seated position.}$ 

# ZG Power Healthcare Recliner | Patient Transfer

# **Select Aesthetic Style** (Denotes back & arm design)










Single Fold Down Arm (Koncert style pictured)

# **Functional Style** (Aesthetic style code precedes selected model number below, e.g., KCZP17)

#### CZP67 ZG Power Healthcare Recliner - Single Fold Down Arm

СОМ	Dillon/Carson	$A_5$	B <sub>6</sub>	$C_7$	D <sub>8</sub>	E9	F <sub>10</sub>	G <sub>11</sub>	H <sub>12</sub>	I <sub>13</sub>
10129	10399	10659	10919	11179	11449	11839	12239	12759	13299	13949

#### **Standard Features**







Recline Functions Operated by Intuitive Control Wand

See pg 2 for finishes



4" Central Locking Caster System



Weightless Zero Gravity Position via Wand Preset

PATIENT WELLNESS & COMFORT OPTIONS

NOTE: Transfer arm standard on right side (facing). Left side (facing) available, please specify if requested.

# **Dimensions**

	Overall —			Seat —			Kick Out				Distance
	W	D	Н	W	D	Н	Back H	Arm H	Seat D	Recline D	from Wall
RCZP67   KCZP67	33	37	46.5	24	21	22	27.5	26.5	40.5	68	15.5
KMCZP67	33	37	46.5	24	21	22	27.5	26.5	40.5	68	15.5

General Info

COM: 7 yds
Static
Wt Cap: 650 lb
Dynamic
Wt Cap: 500 lb
Unit
Weight: 174 lb

## **Available Options** \*Refer to pg 99 for full option details and images\*

ARM CAP OPTIONS

	THAT CIT OF FIGURE	2 for finishes		TATIETT WELLINESS & CONTROLL OF THOS	15
U	Urethane	Standard	ABP	Adjustable Back Pad	
J1	Solid Surface - Group 1	619		Textile Grades COM-B	459
J2	Solid Surface - Group 2	819		Textile Grades C-I	569
			AHP	Adjustable Head Rest Pillow	339
			DP	Dual IV Poles & Brackets <sup>1</sup>	498
			FBH	Foley Bag Holder	89
	CASTER OPTIONS		FT	Foot Tray	279
Z-Q3	3" - Dual Wheel Central Locking System	n No Charge	IV	IV Pole & Bracket <sup>1</sup>	249
Z-Q3	4" - Central Locking System	Standard	LC	LiquiCell	229
Z-Q4 Z-Q5	5" - Central Locking System	119	LCKFR	Orthopedic Locking Footrest <sup>6</sup>	Standard
Z-Q3	5 - Central Locking System	117	MBB	Moisture Barrier (Back)	59
			MBS	Moisture Barrier (Seat)	59
			MR	Monitor Arm (right side facing) <sup>1,2</sup>	849
	HEAT/MASSAGE OPTIONS (IEC 6060)	1_1 CEDTIEIED)	MTR	Monitor Arm with Tray (right side facing) <sup>1,2</sup>	949
	TIEAT/ MASSAGE OF HONS (IEC 0000	I-I CERTIFIED)	PR1	Paper Roll Holder <sup>3</sup>	229
	Wand Controls		TE	Trendelenburg	Standard
HW	Heat	1139	USB	USB Power Outlet (IEC 60601-1 Certified)	459
HMW	Heat & Massage	1369	WR	Wrist Restraint Mounts (pair)	229
	In-Arm Controls			TRAY TABLE OPTIONS	
Н	Heat	849		TRAT TABLE OF TIONS	
HM	Heat & Massage	1119	TSL1	Adjustable (Left) Side Tray, on Stationary A	
			TSR167	Adjustable (Right) Side Tray, on Fold Down A	\rm <sup>1,2,5</sup> 449
	ADDITIONAL FEATURES & OPTIONS		LFSTR	Large Folding (Right) Side Tray 2,4	249
BBU	Battery Back Up	649	LFSTRCH	Large Folding Cup Holder (Right) Side Tray	319
PB	Push Bar	Standard	PLFSTR	Padded Large Folding (Right) Side Tray <sup>2</sup>	379
EG	Protective Edge Guard	Standard	SFSTR	Small Folding (Right) Side Tray 2,4	249
QR CAL	Additional Quick Release Accessory Bra CAL TB133			FOLD DOWN TRANSFER ARM SPECIFIC	
Z	Multi Cover Use Pricing for	Highest Grade	НС	Hinge Cover	Standard

<sup>&</sup>lt;sup>1</sup>Includes QR bracket connection installed on recliner

<sup>&</sup>lt;sup>6</sup>Position is locked in place by motors and can only be adjusted by the control wand.







<sup>&</sup>lt;sup>2</sup>Standard on right side (facing). Replace "R" with "L" for left side (facing) option where available.

 $<sup>^3</sup> Takes\ place\ of\ Push\ Bar.\ Paper\ sold\ separately.$ 

<sup>&</sup>lt;sup>4</sup>Contact customer service for solid surface pricing

 $<sup>{}^{\</sup>scriptscriptstyle 5}\!\text{Over the lap tray function recommended in seated position.}$ 

# ZG Power Healthcare Recliner | Patient Transfer

# **Select Aesthetic Style** (Denotes back & arm design)









Dual Fold Down Arms (Koncert style pictured)

# **Functional Style** (Aesthetic style code precedes selected model number below, e.g., KCZP17)

# CZP62 ZG Power Healthcare Recliner - Dual Fold Down Arms

СОМ	Dillon/Carson	$A_5$	B <sub>6</sub>	$C_7$	D <sub>8</sub>	E9	F <sub>10</sub>	G <sub>11</sub>	$H_{12}$	I <sub>13</sub>
10789	11049	11319	11579	11839	12109	12499	12899	13419	13949	14599

#### **Standard Features**







Recline Functions Operated by Intuitive Control Wand



4" Central Locking Caster System



Weightless Zero Gravity Position via Wand Preset

Dimensions			──Overall ── W D H			W D H			Arm H	Kick Out Seat D Recline D		Distance from Wall	
	RCZP62   KCZP62	33	37	46.5	24	21	22	27.5	26.5	40.5	68	15.5	
	KMCZP62	33	37	46.5	24	21	22	27.5	26.5	40.5	68	15.5	

General Info
COM: 7 yds
Static
Wt Cap: 650 lb
Dynamic
Wt Cap: 500 lb
Unit
Weight: 174 lb

## **Available Options** \*Refer to pg 99 for full option details and images\*

	ARM CAP OPTIONS See pg 2	for finishes		PATIENT WELLNESS & COMFORT OPTION	IS
U	Urethane	Standard	ABP	Adjustable Back Pad	
J1	Solid Surface - Group 1	619		Textile Grades COM-B	459
J2	Solid Surface - Group 2	819		Textile Grades C-I	569
			AHP	Adjustable Head Rest Pillow	339
			DP	Dual IV Poles & Brackets <sup>1</sup>	498
			FBH	Foley Bag Holder	89
	CASTER OPTIONS		FT	Foot Tray	279
Z-Q3	3" Dual Wheel Central Leaking System	No Chargo	IV	IV Pole & Bracket <sup>1</sup>	249
Z-Q3 Z-Q4	3" - Dual Wheel Central Locking System 4" - Central Locking System	No Charge Standard	LC	LiquiCell	229
Z-Q4 Z-Q5	5" - Central Locking System	119	LCKFR	Orthopedic Locking Footrest <sup>6</sup>	Standard
Z-Q5	5 - Central Locking System	119	MBB	Moisture Barrier (Back)	59
			MBS	Moisture Barrier (Seat)	59
			MR	Monitor Arm (right side facing) <sup>1,2</sup>	849
	HEAT/MASSAGE OPTIONS (IEC 60601-	1 CEDTIEIED)	MTR	Monitor Arm with Tray (right side facing) <sup>1,2</sup>	949
	TIEAT/ MASSAGE OF HONS (IEC 60001-	(CERTIFIED)	PR1	Paper Roll Holder <sup>3</sup>	229
	Wand Controls		TE	Trendelenburg	Standard
HW	Heat	1139	USB	USB Power Outlet (IEC 60601-1 Certified)	459
HMW	Heat & Massage	1369	WR	Wrist Restraint Mounts (pair)	229
				TRAY TABLE OPTIONS	
			TD362	Dual Adjustable Side Trays, on Fold Down A	rms <sup>1,5</sup> 898
			TSR167	Adjustable (Right) Side Tray, on Fold Down A	rm <sup>1,5</sup> 449
	ADDITIONAL FEATURES & OPTIONS		TSL167	Adjustable (Left) Side Tray, on Fold Down A	rm <sup>1,5</sup> 449
BBU	Battery Back Up	649			
PB	Push Bar	Standard			
EG	Protective Edge Guard	Standard			
QR	Additional Quick Release Accessory Brac	ket 69		FOLD DOWN TRANSFER ARM SPECIFIC	
CAL	CAL TB133	See pg 9		FOLD DOWN TRAINSFER ARIVI SPECIFIC	
Z	Multi Cover Use Pricing for H	lighest Grade	HC	Hinge Cover	Standard

<sup>&</sup>lt;sup>1</sup>Includes QR bracket connection installed on recliner

<sup>&</sup>lt;sup>6</sup>Position is locked in place by motors and can only be adjusted by the control wand.







<sup>&</sup>lt;sup>2</sup>Standard on right side (facing). Replace "R" with "L" for left side (facing) option where available.

 $<sup>^3</sup> Takes\ place\ of\ Push\ Bar.\ Paper\ sold\ separately.$ 

<sup>&</sup>lt;sup>4</sup>Contact customer service for solid surface pricing

<sup>&</sup>lt;sup>5</sup>Over the lap tray function recommended in seated position.

# ZG Lift Assist Power Healthcare Recliner

# **Select Aesthetic Style** (Denotes back & arm design)









(Koncert style pictured)

**Functional Style** (Aesthetic style code precedes selected model number below, e.g., KCZL17)

#### CZL17 ZG Lift Assist Power Healthcare Recliner

COM	Dillon/Carson	$A_5$	В <sub>6</sub>	$C_7$	$D_8$	E9	F <sub>10</sub>	G <sub>11</sub>	H <sub>12</sub>	I <sub>13</sub>
11199	11999	12499	12899	13029	13159	13419	13699	14349	14749	15269

#### **Standard Features**







Recline Functions Operated by Intuitive Control Wand



3" Automatic Central Locking Caster System



Weightless Zero Gravity Position via Wand Preset



Lift Assist Function Aids User in Exiting Chair

# **Dimensions**

	Overall —				– Seat -		Kick Out				Distance
	W	D	Н	W	D	Н	Back H	Arm H	Seat D	Recline D	from Wall
RCZL17   KCZL17	30.5	38.5	47.5	24	21.5	23.5	28	29.5	39	87.5	18.75
KMCZL17	30.5	38.5	47.5	24	21.5	23.5	28	29.5	39	87.5	18.75



# **Available Options**

ARM CAP OPTIONS	See pg 2 for finishes
Urethane	169

 U
 Urethane
 169

 J1
 Solid Surface - Group 1
 619

 J2
 Solid Surface - Group 2
 819

 F
 Fully Upholstered Arms
 Standard

#### **CASTER OPTIONS**

ZL-Q3	3" Auto-Locking Central Caster System	Standard
ZL-Q4	4" Auto-Locking Central Caster System Increases seat height by 1"	No Charge
ZL-Q5	5" Auto-Locking Central Caster System	119

## HEAT/MASSAGE OPTIONS (IEC 60601-1 CERTIFIED)

	Wand Controls	
HW	Heat	1139
HMW	Heat & Massage	1369
	In-Arm Controls	
Н	Heat	849
НМ	Heat & Massage	1119

#### ADDITIONAL FEATURES & OPTIONS

BBU	Battery Back Up	649
PB	Push Bar	Standard
EG	Protective Edge Guard	Standard
CAL	CAL TB133	See pg 9
Z	Multi Cover	Use Pricing for Highest Grade

# PATIENT WELLNESS & COMFORT OPTIONS

ABP	Adjustable Back Pad	
	Textile Grades COM-B	459
	Textile Grades C-I	569
AHP	Adjustable Head Rest Pillow	339
FT	Foot Tray	279
LC	LiquiCell	229
LCKFR	Orthopedic Locking Footrest <sup>6</sup>	Standard
MBB	Moisture Barrier (Back)	59
MBS	Moisture Barrier (Seat)	59
PR1	Paper Roll Holder <sup>3</sup>	229
TE	Trendelenburg	Standard
USB	USB Power Outlet (IEC 60601-1 Certified)	459







<sup>\*</sup>Refer to pg 99 for full option details and images\*

<sup>&</sup>lt;sup>3</sup>Takes place of Push Bar. Paper sold separately.

<sup>&</sup>lt;sup>6</sup>Position is locked in place by motors and can only be adjusted by the control wand.

# **DURABLE® POWER SERIES** ZG Lift Assist Power Healthcare Recliner | Patient Transfer

# **Select Aesthetic Style** (Denotes back & arm design)







**Functional Style** (Aesthetic style code precedes selected model number below, e.g., KCZL17)

#### CZL67 ZG Lift Assist Power Healthcare Recliner - Single Fold Down Arm

COM Dillon/Carson  $C_7$ D<sub>8</sub> E9 F<sub>10</sub>  $G_{11}$  $H_{12}$ I<sub>13</sub> 12899 13159 13429 13699 13999 14479 14869 15399 15999



Single Fold Down Arm (Koncert style pictured)

**NOTE:** Transfer arm standard on right side (facing). Left side (facing) available, please specify if requested.

#### **Standard Features**





ARM CAP OPTIONS

In-Arm Controls

Heat & Massage

Battery Back Up

Protective Edge Guard

Push Bar

CAL TB133

Multi Cover

H Heat

НМ

BBU

PB

CAL

Ζ



Recline Functions Operated by Intuitive Control Wand

See pg 2 for finishes



3" Automatic Central Locking Caster System



Weightless Zero Gravity Position via Wand Preset

PATIENT WELLNESS & COMFORT OPTIONS



Lift Assist Function Aids User in Exiting Chair

# **Dimensions**

	Overall —		Seat —				Kick Out			Distance	
	W	D	Н	W	D	Н	Back H	Arm H	Seat D	Recline D	from Wall
RCZL67   KCZL67	30.5	38.5	47.5	24	21.5	23.5	28	29.5	39	87.5	18.75
KMCZL67	30.5	38.5	47.5	24	21.5	23.5	28	29.5	39	87.5	18.75

# General Info COM: 7 yds Static Wt Cap: 650 lb Dynamic Wt Cap: 500 lb Unit Weight: 214 lb

229 Standard

459

## **Available Options**

	, ,			
U J1	Urethane Solid Surface - Group 1	Standard 619	ABP	Adjustable Back Pad Textile Grades COM-B
J2	Solid Surface - Group 2	819	AHP FT LC	Textile Grades C-I Adjustable Head Rest Pillow Foot Tray LiquiCell
	CASTER OPTIONS		LCKFR	Orthopedic Locking Footrest <sup>6</sup>
ZL-Q3 ZL-Q4 ZL-Q5	3" Auto-Locking Central Caster System 4" Auto-Locking Central Caster System Increases seat height by 1"  5" Auto-Locking Central Caster System Increases seat height by 2"	Standard No Charge 119	MBB MBS PR1 TE USB	Moisture Barrier (Back) Moisture Barrier (Seat) Paper Roll Holder <sup>3</sup> Trendelenburg USB Power Outlet (IEC 60601-1 Certified)
	HEAT/MASSAGE OPTIONS (IEC 60601-1	CERTIFIED)		
HW HMW	Wand Controls Heat Heat & Massage	1139 1369		

849

1119

649

Standard

Standard

See pg 9

Use Pricing for Highest Grade

**ADDITIONAL FEATURES & OPTIONS** 







<sup>\*</sup>Refer to pg 99 for full option details and images\*

<sup>&</sup>lt;sup>3</sup>Takes place of Push Bar. Paper sold separately.

<sup>&</sup>lt;sup>6</sup>Position is locked in place by motors and can only be adjusted by the control wand.

# **DURABLE® POWER SERIES** ZG Lift Assist Power Healthcare Recliner | Patient Transfer

# **Select Aesthetic Style** (Denotes back & arm design)










Dual Fold Down Arms (Koncert style pictured)

# **Functional Style** (Aesthetic style code precedes selected model number below, e.g., KCZL17)

#### CZL62 ZG Lift Assist Power Healthcare Recliner - Dual Fold Down Arms

СОМ	Dillon/Carson	$A_5$	B <sub>6</sub>	$C_7$	D <sub>8</sub>	E9	F <sub>10</sub>	G <sub>11</sub>	H <sub>12</sub>	I <sub>13</sub>
12869	13159	13299	13679	13879	14209	14479	14869	15399	15799	16319

#### **Standard Features**







Recline Functions Operated by Intuitive Control Wand



3" Automatic Central Locking Caster System



Weightless Zero Gravity Position via Wand Preset



Lift Assist Function Aids User in Exiting Chair

# **Dimensions**

		- Overall -		Seat —		Kick Out			Distance		
	W	D	Н	W	D	Н	Back H	Arm H	Seat D	Recline D	from Wall
RCZL62   KCZL62	30.5	38.5	47.5	24	21.5	23.5	28	29.5	39	87.5	18.75
KMCZL62	30.5	38.5	47.5	24	21.5	23.5	28	29.5	39	87.5	18.75

Genera	ıl Info
сом:	7 yds
Static Wt Cap:	650 lb
Dynamic <b>Wt Cap:</b>	500 lb
Unit <b>Weight:</b>	214 lb

# **Available Options**

	ARM CAP OPTIONS	See pg 2 for finishes
U	Urethane	Standard
J1	Solid Surface - Group 1	619
J2	Solid Surface - Group 2	819

#### **CASTER OPTIONS**

ZL-Q3	3" Auto-Locking Central Caster System	Standard
ZL-Q4	4" Auto-Locking Central Caster System Increases seat height by 1"	No Charge
ZL-Q5	5" Auto-Locking Central Caster System	119

## HEAT/MASSAGE OPTIONS (IEC 60601-1 CERTIFIED)

	Wand Controls	
HW	Heat	1139
HMW	Heat & Massage	1369

#### ADDITIONAL FEATURES & OPTIONS

BBU	Battery Back Up	649
PB	Push Bar	Standard
EG	Protective Edge Guard	Standard
CAL	CAL TB133	See pg 9
Z	Multi Cover	Use Pricing for Highest Grade

# PATIENT WELLNESS & COMFORT OPTIONS

ABP	Adjustable Back Pad	
	Textile Grades COM-B	459
	Textile Grades C-I	569
AHP	Adjustable Head Rest Pillow	339
FT	Foot Tray	279
LC	LiquiCell	229
LCKFR	Orthopedic Locking Footrest <sup>6</sup>	Standard
MBB	Moisture Barrier (Back)	59
MBS	Moisture Barrier (Seat)	59
PR1	Paper Roll Holder <sup>3</sup>	229
TE	Trendelenburg	Standard
USB	USB Power Outlet (IEC 60601-1 Certified)	459







<sup>\*</sup>Refer to pg 99 for full option details and images\*

<sup>&</sup>lt;sup>3</sup>Takes place of Push Bar. Paper sold separately.

<sup>&</sup>lt;sup>6</sup>Position is locked in place by motors and can only be adjusted by the control wand.

# ZG Exam Lift Power Healthcare Recliner

# **Select Aesthetic Style** (Denotes back & arm design)







(Koncert style pictured)

**Functional Style** (Aesthetic style code precedes selected model number below, e.g., KCZX67)

#### CZX17 ZG Exam Lift Power Healthcare Recliner

СОМ	Dillon/Carson	$A_5$	B <sub>6</sub>	$C_7$	D <sub>8</sub>	E9	F <sub>10</sub>	G <sub>11</sub>	H <sub>12</sub>	I <sub>13</sub>
16319	16579	16709	17109	17369	17629	17899	18419	18949	19469	19869

#### **Standard Features**





ARM CAP OPTIONS



Recline Functions Operated by Intuitive Control Wand

See pg 2 for finishes



3" Dual Wheel Central Locking Caster System



Weightless Zero Gravity Position via Wand Preset

PATIENT WELLNESS & COMFORT OPTIONS



Raise to Exam Height via Non-Handed Foot Pedals

# **Dimensions**

	Overall				Seat —			Kick Out			
	W	D	Н	W	D	Н	Back H	Arm H	Seat D	Recline D	from Wall
RCZX17   KCZX17	34.5	38	48	24	21	23	26	28	40	73	20
KMCZX17	34.5	38	48	24	21	23	26	28	40	73	20

General Info				
сом:	7 yds			
Static Wt Cap:	650 lb			
Dynamic <b>Wt Cap:</b>	500 lb			
Unit <b>Weight:</b>	281 lb			

## **Available Options**

	711111 C/11 O/ 110113 Sec pg 2 /				_
U	Urethane	169	ABP	Adjustable Back Pad	
J1	Solid Surface - Group 1	619		Textile Grades COM-B	459
J2	Solid Surface - Group 2	819		Textile Grades C-I	569
F	Fully Upholstered Arms	Standard	AHP	Adjustable Head Rest Pillow	339
			DP	Dual IV Poles & Brackets <sup>1</sup>	498
			FBH	Foley Bag Holder	89
	CASTER OPTIONS		IV	IV Pole & Bracket <sup>1</sup>	249
Z-Q3	3" - Dual Wheel Central Locking Casters	Standard	LC	LiquiCell	229
Z-Q3 Z-Q4	4" - Central Locking System	No Charge	LCKFR	Orthopedic Locking Footrest <sup>6</sup>	Standard
Z-Q4	Increases seat height by 1".	No Charge	MBB	Moisture Barrier (Back)	59
Z-Q5	5" - Central Locking System	119	MBS	Moisture Barrier (Seat)	59
<u> 2</u> -Q3	Increases seat height by 2".	117	MR	Monitor Arm (right side facing) <sup>1,2</sup>	849
	mercases seat neight by 2 .		MTR	Monitor Arm with Tray (right side facing) <sup>1,2</sup>	949
	HEAT/MASSAGE OPTIONS (IEC 60601-1	CERTIFIED)	PR1	Paper Roll Holder <sup>3</sup>	229
	TIEMI/ MIASSAGE OF TIONS (IEC 00001-1	CERTIFIED)	TE	Trendelenburg	Standard
	Wand Controls		USB	USB Power Outlet (IEC 60601-1 Certified)	459
HW	Heat	1139	WR	Wrist Restraint Mounts (pair)	229
HMW	Heat & Massage	1369			
	In-Arm Controls				
Н	Heat	849		TRAY TABLE OPTIONS	
НМ	Heat & Massage	1119	TSR1	Adjustable (Right) Side Tray <sup>1,2,5</sup>	399
	9		TD3	Dual Adjustable Side Trays (TSR1 & TSL1) <sup>1,5</sup>	798
	ADDITIONAL FEATURES & OPTIONS		LFSTR	Large Folding (Right) Side Tray <sup>4</sup>	249
DDLI	D. II. D. I. II.	(40	LFSTRCH	Large Folding Cup Holder (Right) Side Tray	319
BBU	Battery Back Up	649	PLFSTR	Padded Large Folding (Right) Side Tray	379
PB	Push Bar	Standard	SFSTR	Small Folding (Right) Side Tray <sup>4</sup>	249
EG	Protective Edge Guard	Standard			
QR	Additional Quick Release Accessory Brack		+0		
CAL	CAL TB133	See pg 9	^Re	efer to pg 99 for full option details and images*	
Z	Multi Cover Use Pricing for H	ignest Grade			

<sup>&</sup>lt;sup>1</sup>Includes QR bracket connection installed on recliner

<sup>&</sup>lt;sup>6</sup>Position is locked in place by motors and can only be adjusted by the control wand.







 $<sup>^2</sup>Standard\ on\ right\ side\ (facing).$  Replace "R" with "L" for left side (facing) option where available.

 $<sup>^3</sup> Takes\ place\ of\ Push\ Bar.\ Paper\ sold\ separately.$ 

<sup>&</sup>lt;sup>4</sup>Contact customer service for solid surface pricing

 $<sup>^{\</sup>rm 5}{\rm Over}$  the lap tray function recommended in seated position.

# **DURABLE® POWER SERIES** ZG Exam Lift Power Healthcare Recliner | Patient Transfer

# **Select Aesthetic Style** (Denotes back & arm design)







Functional Style (Aesthetic style code precedes selected model number below, e.g., KCZX67)

#### CZX67 ZG Exam Lift Power Healthcare Recliner - Single Fold Down Arm

COM Dillon/Carson A<sub>5</sub>  $C_7$ E9 F<sub>10</sub>  $G_{11}$  $H_{12}$ I<sub>13</sub> 17109 17239 17629 17899 18159 18549 19079 19599 20129 20659



Single Fold Down Arm (Koncert style pictured)

**NOTE:** Transfer arm standard on right side (facing). Left side (facing) available, please specify if requested.

#### **Standard Features**







Recline Functions Operated by Intuitive Control Wand



3" Dual Wheel Central Locking Caster System



Weightless Zero Gravity Position via Wand Preset



Raise to Exam Height via Non-Handed Foot Pedals

## **Dimensions**

		Overall —		Seat -			Kick Out			Distance		
		W	D	Н	W	D	Н	Back H	Arm H	Seat D	Recline D	from Wall
RCZX67   KG	CZX67	34.5	38	48	24	21	23	26	28	40	73	20
KM	CZX67	34.5	38	48	24	21	23	26	28	40	73	20

General Info

COM: 8 yds
Static
Wt Cap: 650 lb
Dynamic
Wt Cap: 500 lb
Unit
Weight: 281 lb

## **Available Options** \*Refer to pg 99 for full option details and images\*

ARM CAP OPTIONS	See pg 2 for finishes

		11.75
U	Urethane	Standard
J1	Solid Surface - Group 1	619
J2	Solid Surface - Group 2	819

#### **CASTER OPTIONS**

Z-Q3	3" - Dual Wheel Central Locking Casters	Standard
Z-Q4	4" - Central Locking System Increases seat height by 1".	No Charge
Z-Q5	5" - Central Locking System Increases seat height by 2".	119

## HEAT/MASSAGE OPTIONS (IEC 60601-1 CERTIFIED)

HW HMW	Wand Controls Heat Heat & Massage	1139 1369
	In-Arm Controls	
Н	Heat	849
HM	Heat & Massage	1119

#### ADDITIONAL FEATURES & OPTIONS

BBU	Battery Back Up		649
PB	Push Bar		Standard
EG	Protective Edge Guard		Standard
QR	Additional Quick Release	Accessory Bracket	69
CAL	CAL TB133		See pg 9
Z	Multi Cover	Use Pricing for High	est Grade

# PATIENT WELLNESS & COMFORT OPTIONS

Adjustable Back Pad	
Textile Grades COM-B	459
Textile Grades C-I	569
Adjustable Head Rest Pillow	339
Dual IV Poles & Brackets <sup>1</sup>	498
Foley Bag Holder	89
IV Pole & Bracket <sup>1</sup>	249
LiquiCell	229
Orthopedic Locking Footrest <sup>6</sup>	Standard
Moisture Barrier (Back)	59
Moisture Barrier (Seat)	59
Monitor Arm (right side facing) <sup>1,2</sup>	849
Monitor Arm with Tray (right side facing) <sup>1,2</sup>	949
Paper Roll Holder <sup>3</sup>	229
Trendelenburg	Standard
USB Power Outlet (IEC 60601-1 Certified)	459
Wrist Restraint Mounts (pair)	229
	Textile Grades COM-B Textile Grades C-I Adjustable Head Rest Pillow Dual IV Poles & Brackets¹ Foley Bag Holder IV Pole & Bracket¹ LiquiCell Orthopedic Locking Footrest6 Moisture Barrier (Back) Moisture Barrier (Seat) Monitor Arm (right side facing)¹² Monitor Arm with Tray (right side facing)¹² Paper Roll Holder³ Trendelenburg USB Power Outlet (IEC 60601-1 Certified)

## TRAY TABLE OPTIONS

TSL1	Adjustable (Left) Side Tray, on Stationary Arm <sup>1,5</sup>	399
TSR167	Adjustable (Right) Side Tray, on Fold Down Arm <sup>1,2,5</sup>	449
LFSTR	Large Folding (Right) Side Tray 2,4	249
LFSTRCH	Large Folding Cup Holder (Right) Side Tray <sup>2</sup>	319
PLFSTR	Padded Large Folding (Right) Side Tray <sup>2</sup>	379
SFSTR	Small Folding (Right) Side Tray 2,4	249

# FOLD DOWN TRANSFER ARM SPECIFIC

<sup>4</sup>Contact customer service for solid surface pricing <sup>5</sup>Over the lap tray function recommended in seated position.



HC Hinge Cover





Standard

<sup>&</sup>lt;sup>6</sup>Position is locked in place by motors and can only be adjusted by the control wand.

¹Includes QR bracket connection installed on recliner

<sup>&</sup>lt;sup>2</sup>Standard on right side (facing). Replace "R" with "L" for left side (facing) option where available.

 $<sup>^{3} \</sup>mbox{\it Takes}$  place of Push Bar. Paper sold separately.

# **DURABLE® POWER SERIES** ZG Exam Lift Power Healthcare Recliner | Patient Transfer

# **Select Aesthetic Style** (Denotes back & arm design)









Dual Fold Down Arms (Koncert style pictured)

# Functional Style (Aesthetic style code precedes selected model number below, e.g., KCZX67)

#### CZX62 ZG Exam Lift Power Healthcare Recliner - Dual Fold Down Arms

COM Dillon/Carson  $C_7$ E9 F<sub>10</sub>  $G_{11}$  $H_{12}$ I<sub>13</sub> 17629 18159 18419 18679 19079 19739 20259 20919 21579

#### **Standard Features**







Recline Functions Operated by Intuitive Control Wand



3" Dual Wheel Central Locking Caster System



Weightless Zero Gravity Position via Wand Preset



Raise to Exam Height via Non-Handed Foot Pedals

# **Dimensions**

		Overall -			- Seat -				Kick Out		Distance
	W	D	Н	W	D	Н	Back H	Arm H	Seat D	Recline D	from Wall
RCZX62   KCZX62	34.5	38	48	24	21	23	26	28	40	73	20
KMCZX62	34.5	38	48	24	21	23	26	28	40	73	20

Genera	ıl Info
сом:	8 yds
Static Wt Cap:	650 lb
Dynamic <b>Wt Cap:</b>	500 lb
Unit Weight:	281 lb

# Available Options

\*Refer to pg 99 for full option details and images\*

	ARM CAP OPTIONS	See pg 2 for finishes
U J1 J2	Urethane Solid Surface - Group 1 Solid Surface - Group 2	Standard 619 819
	CASTER OPTIONS	
Z-Q3 Z-Q4	3" - Dual Wheel Central Locking 4" - Central Locking System Increases seat height by 1".	g Casters Standard No Charge
Z-Q5	5" - Central Locking System Increases seat height by 2".	119
	HEAT/MASSAGE OPTIONS (IE	C 60601-1 CERTIFIED)

# ADDITIONAL FEATURES & OPTIONS

Wand Controls

Heat & Massage

HW

HMW

BBU	Battery Back Up		649
PB	Push Bar		Standard
EG	Protective Edge Guard		Standard
QR	Additional Quick Releas	e Accessory Bracket	69
CAL	CAL TB133		See pg 9
Z	Multi Cover	Use Pricing for High	est Grade

# PATIENT WELLNESS & COMFORT OPTIONS

ABP	Adjustable Back Pad	
	Textile Grades COM-B	459
	Textile Grades C-I	569
AHP	Adjustable Head Rest Pillow	339
DP	Dual IV Poles & Brackets <sup>1</sup>	498
FBH	Foley Bag Holder	89
IV	IV Pole & Bracket <sup>1</sup>	249
LC	LiquiCell	229
LCKFR	Orthopedic Locking Footrest <sup>6</sup>	Standard
MBB	Moisture Barrier (Back)	59
MBS	Moisture Barrier (Seat)	59
MR	Monitor Arm (right side facing) <sup>1,2</sup>	849
MTR	Monitor Arm with Tray (right side facing) <sup>1,2</sup>	949
PR1	Paper Roll Holder <sup>3</sup>	229
TE	Trendelenburg	Standard
USB	USB Power Outlet (IEC 60601-1 Certified)	459
WR	Wrist Restraint Mounts (pair)	229

## TRAY TABLE OPTIONS

10362	Duai Adjustable Side Trays, on Fold Down Arms	898
TSR167	Adjustable (Right) Side Tray, on Fold Down Arm <sup>1,5</sup>	449
TSL167	Adjustable (Left) Side Tray, on Fold Down Arm <sup>1,5</sup>	449

#### FOLD DOWN TRANSFER ARM SPECIFIC

HC Hinge Cover Standard

<sup>&</sup>lt;sup>6</sup>Position is locked in place by motors and can only be adjusted by the control wand.







1139

1369

<sup>&</sup>lt;sup>1</sup>Includes QR bracket connection installed on recliner

<sup>&</sup>lt;sup>2</sup>Standard on right side (facing). Replace "R" with "L" for left side (facing) option where available.

 $<sup>^3</sup>$ Takes place of Push Bar. Paper sold separately.

<sup>&</sup>lt;sup>4</sup>Contact customer service for solid surface pricing <sup>5</sup>Over the lap tray function recommended in seated position. <sup>6</sup>Position is locked in place by motors and can only be adjuste