



EV, Rêma | Patient Transfer Arm Mobile Rocker Recliners - 24" Seat Width - 550 lb Weight Rating



EV | E1127R | EC1127R
32.5" Overall Width



Rêma | R1127R | RC1127R
32.5" Overall Width

DIMENSIONS | E1127R | Chaise Models

Width	32.5	Seat Width	23.75	Back Height	25
Depth	38	Seat Depth	21.5	Arm Height	25
Height	44.75	Seat Height	19.75	Weight	151
COM (Standard)	6 YDS	COM (Chaise)	7 YDS	Seat Depth (kick-out)	37.5
Reclined Depth	72.5	Weight Capacity	550 lbs.	Distance from Wall	20.5

DIMENSIONS | R1127R | Chaise Models

Width	32.5	Seat Width	23.75	Back Height	23.75
Depth	38	Seat Depth	21.5	Arm Height	25
Height	43.5	Seat Height	19.75	Weight	153
COM (Standard)	6 YDS	COM (Chaise)	7 YDS	Seat Depth (kick-out)	37.5
Reclined Depth	72.5	Weight Capacity	550 lbs.	Distance from Wall	19.5

1127R MOBILE ROCKER RECLINER INFINITE POSITION LOCKING BACK - RIGHT FACING FOLD DOWN ARM - URETHANE ARM CAPS

COM	4	5	6	7	8	9	10	11	
	4839	4899	4909	5039	5209	5289	5569	5819	6169

C1127R MOBILE CHAISE ROCKER RECLINER INFINITE POSITION LOCKING BACK - RIGHT FACING FOLD DOWN ARM - URETHANE ARM CAPS

COM	4	5	6	7	8	9	10	11	
	5139	5279	5309	5379	5589	5739	5979	6279	6599

*If patient transfer arm is specified on left side (facing), then handle will be placed on right side (facing) and vice versa.

Options

U	Urethane Arm Caps*
W	Wood Arm Caps*
J	Solid Surface Arm Caps* - Group 1
Z	Multi Cover


List Price

Standard
\$229
\$619 (Group 2 add \$200)
Use Highest Grade


Options

CAL	CAL TB133	See page 12
MBB	Moisture Barrier Back	\$59
MBS	Moisture Barrier Seat	\$59
FBH	Foley Bag Holder	\$89
3S4	All 3" Swivel Locking Caster	\$19
AHP	Adjustable Headrest Pillow (*Requires 1 yd of fabric)	\$339
QR	Additional Quick Release Accessory Bracket	\$69
MR/L	Right/Left Side (facing) Monitor Arm (non-fold down arm only)	\$629
MTR/L	Right/Left Side (facing) Monitor Arm with Tray (non-fold down arm only)	\$629
ABP	Adjustable Back Pad (COM - Grade 6)	\$459
	Adjustable Back Pad (Grade 7-11)	\$569

*Please refer to page 14 for finishes.



Infinite Position Locking Back Lever
(1127 + 4127 Models Only)
Infinite Position Locking Back lever allows patient/caregivers to lock back in any recline position. Locking Back Lever is located on the inside of the left arm (facing) on select models. Locking Back Lever is located on the outside of the left arm (facing) on Dual Lateral Transfer Arm models.



Rocker/Lock Mechanism
(1127 + 4127 Models Only)
The mobile rocker recliner features a locking mechanism that locks the rocker making it easy on the patient for entering/exiting the chair and also during transport.

Please refer to pages 134/135 for Option descriptions & images.
Please refer to page 132 for the recliner ordering guide.

