



EV, Rêma | Patient Transfer Arm Mobile Rocker Recliners - 24" Seat Width - 550 lb Weight Rating



EV | E1127L | EC1127L
32.5" Overall Width



Rêma | R1127L | RC1127L
32.5" Overall Width

DIMENSIONS | E1127L | Chaise Models

Width	32.5	Seat Width	23.75	Back Height	25
Depth	38	Seat Depth	21.5	Arm Height	25
Height	44.75	Seat Height	19.75	Weight	151
COM (Standard)	6 YDS	COM (Chaise)	7 YDS	Seat Depth (kick-out)	37.5
Reclined Depth	72.5	Weight Capacity	550 lbs.	Distance from Wall	20.5

DIMENSIONS | R1127L | Chaise Models

Width	32.5	Seat Width	23.75	Back Height	23.75
Depth	38	Seat Depth	21.5	Arm Height	25
Height	43.5	Seat Height	19.75	Weight	153
COM (Standard)	6 YDS	COM (Chaise)	7 YDS	Seat Depth (kick-out)	37.5
Reclined Depth	72.5	Weight Capacity	550 lbs.	Distance from Wall	19.5

1127L MOBILE ROCKER RECLINER INFINITE POSITION LOCKING BACK - LEFT FACING FOLD DOWN ARM - URETHANE ARM CAPS

COM	4	5	6	7	8	9	10	11
4419	4479	4489	4599	4759	4829	5089	5319	5629


C1127L MOBILE CHAISE ROCKER RECLINER INFINITE POSITION LOCKING BACK - LEFT FACING FOLD DOWN ARM - URETHANE ARM CAPS

COM	4	5	6	7	8	9	10	11
4699	4819	4849	4919	5109	5239	5459	5739	6029

*If patient transfer arm is specified on left side (facing), then handle will be placed on right side (facing) and vice versa.


Options	List Price	Options	List Price
U Urethane Arm Caps*	Standard	MBS Moisture Barrier Seat	\$50
W Wood Arm Caps*	\$199	FBH Foley Bag Holder	\$79
J Solid Surface Arm Caps*	\$539	3S4 All 3" Swivel Locking Caster	\$19
Z Multi Cover	Use Highest Grade	AHP Adjustable Headrest Pillow	\$299
CAL CAL TB133	See page 12	(*Requires 1 yd of fabric)	
MBB Moisture Barrier Back	\$50	QR Quick Release Bracket	\$59
		MR/L Right/Left Side (facing) Monitor Arm	\$549
		(non-fold down arm only)	
		MTR/L Right/Left Side (facing) Monitor Arm with Tray	\$549
		(non-fold down arm only)	

*Please refer to page 14 for finishes.



Infinite Position Locking Back Lever
(1127 + 4127 Models Only)

Infinite Position Locking Back lever allows patient/caregivers to lock back in any recline position. Locking Back Lever is located on the inside of the left arm (facing) on select models. Locking Back Lever is located on the outside of the left arm (facing) on Dual Lateral Transfer Arm models.



Rocker/Lock Mechanism
(1127 + 4127 Models Only)

The mobile rocker recliner features a locking mechanism that locks the rocker making it easy on the patient for entering/exiting the chair and also during transport.

Please refer to pages 150/151 for Options descriptions & images.
Please refer to page 148 for the recliner ordering guide.

