DURABLE® POWER SERIES ZG Exam Lift Power Healthcare Recliner | Patient Transfer

Select Aesthetic Style (Denotes back & arm design)







Functional Style (Aesthetic style code precedes selected model number below, e.g., KCZX67)

CZX67 ZG Exam Lift Power Healthcare Recliner - Single Fold Down Arm

COM Dillon/Carson C_7 E9 F₁₀ G_{11} H₁₂ I₁₃ 17899 18159 18549 16839 17109 17239 17629 19079 19599 20129 20659



Single Fold Down Arm (Koncert style pictured)

NOTE: Transfer arm standard on right side (facing). Left side (facing) available, please specify if requested.

Standard Features



Leg Support

ARM CAP OPTIONS









3" Dual Wheel Central Locking Caster System



Weightless Zero Gravity Position via Wand Preset

PATIENT WELLNESS & COMFORT OPTIONS



Raise to Exam Height via Non-Handed Foot Pedals

	Overall —			Seat —					Kick Out	Distance	
	W	D	Н	W	D	Н	Back H	Arm H	Seat D	Recline D	from Wall
RCZX67 KCZX67	34.5	38	48	24	21	23	26	28	40	73	20
KMCZX67	34.5	38	48	24	21	23	26	28	40	73	20

General Info

Available Options

Dimensions

Refer to pg 103 for full option details and images

See pg 2 for finishes

		/						
U	Urethane	Standard	ABP	Adjustable Back Pad				
J1	Solid Surface - Group 1	619		Textile Grades COM-B	459			
J2	Solid Surface - Group 2	819		Textile Grades C-I	569			
			AHP	Adjustable Head Rest Pillow	339			
			DP	Dual IV Poles & Brackets ¹	498			
			FBH	Foley Bag Holder	89			
	CASTER OPTIONS		IV	IV Pole & Bracket ¹	249			
Z-Q3	3" - Dual Wheel Central Locking Casters	Standard	LC	LiquiCell	229			
Z-Q3 Z-Q4	4" - Central Locking System	No Charge	LCKFR	Orthopedic Locking Footrest ⁶	Standard			
Z-Q4	Increases seat height by 1".	No Charge	MBB	Moisture Barrier (Back)	59			
Z-Q5 5" - Central Loc	, , , , , , , , , , , , , , , , , , ,	119	MBS	Moisture Barrier (Seat)	59			
	Increases seat height by 2".	119	MR	Monitor Arm (right side facing) ^{1,2}	849			
	mercuses seat neight by 2 .		MTR	Monitor Arm with Tray (right side facing) ^{1,2}	949			
	HEAT/MASSAGE OPTIONS (IEC 60601-1	CERTIFIED)	PR1	Paper Roll Holder ³	229			
	HEAT/MASSAGE OF HONS (IEC 60601-	CERTIFIED)	TE	Trendelenburg	Standard			
	Wand Controls		USB	USB Power Outlet (IEC 60601-1 Certified)	459			
HW	Heat	1139	WR	Wrist Restraint Mounts (pair)	229			
HMW	Heat & Massage	1369						
	In-Arm Controls			TRAY TABLE OPTIONS				
Н	Heat	849		TRAT TABLE OF TIONS				
HM	Heat & Massage	1119	TSL1	Adjustable (Left) Side Tray, on Stationary A				
			TSR167	Adjustable (Right) Side Tray, on Fold Down Arm ^{1,2,5} 449				
	ADDITIONAL FEATURES & OPTIONS		LFSTR	Large Folding (Right) Side Tray ^{2, 4}	249			
BBU	Battery Back Up	649	LFSTRCH	Large Folding Cup Holder (Right) Side Tray ²	319			
PB	Push Bar	Standard	PLFSTR	Padded Large Folding (Right) Side Tray ²	379			
EG	Protective Edge Guard	Standard	SFSTR	Small Folding (Right) Side Tray 2,4	249			
QR	Additional Quick Release Accessory Brack							
CAL	CAL TB133	See pg 9		FOLD DOWN TRANSFER ARM SPECIFIC				
Z	Multi Cover Use Pricing for H		нс	Hinge Cover	Standard			
_	Multi Cover Ose Fricing for F	ignest Grade	I IIC	i iiige covei	Standard			

¹Includes QR bracket connection installed on recliner

⁶Position is locked in place by motors and can only be adjusted by the control wand.







²Standard on right side (facing). Replace "R" with "L" for left side (facing) option where available.

³Takes place of Push Bar. Paper sold separately.

⁴Contact customer service for solid surface pricing

⁵Over the lap tray function recommended in seated position.