DURABLE® POWER SERIES

ZG Power Healthcare Recliner | Patient Transfer

Select Aesthetic Style (Denotes back & arm design)










Single Fold Down Arm (Koncert style pictured)

Functional Style (Aesthetic style code precedes selected model number below, e.g., KCZP17)

CZP67 ZG Power Healthcare Recliner - Single Fold Down Arm

СОМ	Dillon/Carson	A_5	B ₆	C_7	D ₈	E9	F ₁₀	G ₁₁	H ₁₂	I ₁₃
10129	10399	10659	10919	11179	11449	11839	12239	12759	13299	13949

Standard Features








Recline Functions Operated by Intuitive Control Wand

See pg 2 for finishes



4" Central Locking Caster System



Weightless Zero Gravity Position via Wand Preset

PATIENT WELLNESS & COMFORT OPTIONS

NOTE: Transfer arm standard on right side (facing). Left side (facing) available, please specify if requested.

Dimensions

	Overall			Seat					Kick Out	Distance	
	W	D	Н	W	D	Н	Back H	Arm H	Seat D	Recline D	from Wall
RCZP67 KCZP67	33	37	46.5	24	21	22	27.5	26.5	40.5	68	15.5
KMCZP67	33	37	46.5	24	21	22	27.5	26.5	40.5	68	15.5

General Info

COM: 7 yds
Static
Wt Cap: 650 lb
Dynamic
Wt Cap: 500 lb
Unit
Weight: 174 lb

Available Options *Refer to pg 103 for full option details and images*

ARM CAP OPTIONS

	THE SEC P	g Z for firmsites		TATILITY WELLINESS & COMITOR OF THOSE				
U	Urethane	Standard	ABP	Adjustable Back Pad				
J1	Solid Surface - Group 1	619		Textile Grades COM-B	459			
J2	Solid Surface - Group 2	819		Textile Grades C-I	569			
			AHP	Adjustable Head Rest Pillow	339			
			DP	Dual IV Poles & Brackets ¹	498			
			FBH	Foley Bag Holder	89			
	CASTER OPTIONS		FT	Foot Tray	279			
7 02	3" Dual Wheel Central Leaking Syste	na Na Chavas	IV	IV Pole & Bracket ¹	249			
Z-Q3	3" - Dual Wheel Central Locking Syste	0	LC	LiquiCell	229			
Z-Q4 Z-Q5	4" - Central Locking System	Standard	LCKFR	Orthopedic Locking Footrest ⁶	Standard			
	5" - Central Locking System	119	MBB	Moisture Barrier (Back)	59			
			MBS	Moisture Barrier (Seat)	59			
			MR	Monitor Arm (right side facing) ^{1,2}	849			
	LIEAT (NAMES ASS OPTIONS (IES COS)	01.1 (EDTIFIED)	MTR	Monitor Arm with Tray (right side facing) ^{1,2}	949			
	HEAT/MASSAGE OPTIONS (IEC 606)	OI-I CERTIFIED)	PR1	Paper Roll Holder ³	229			
	Wand Controls		TE	Trendelenburg	Standard			
HW	Heat	1139	USB	USB Power Outlet (IEC 60601-1 Certified)	459			
HMW	Heat & Massage	1369	WR	Wrist Restraint Mounts (pair)	229			
	In-Arm Controls			TRAY TABLE OPTIONS				
Н	Heat	849		TRAT TABLE OF HOINS				
HM	Heat & Massage	1119	TSL1	Adjustable (Left) Side Tray, on Stationary A	rm ^{1,5} 399			
			TSR167	Adjustable (Right) Side Tray, on Fold Down Arm ^{1,2,5} 449				
	ADDITIONAL FEATURES & OPTIONS		LFSTR	Large Folding (Right) Side Tray ^{2, 4}	249			
BBU	Battery Back Up	649	LFSTRCH	Large Folding Cup Holder (Right) Side Tray	319			
PR	Push Bar	Standard	PLFSTR	Padded Large Folding (Right) Side Tray ²	379			
		Standard	SFSTR	Small Folding (Right) Side Tray 2,4	249			
EG	Protective Edge Guard							
QR	Additional Quick Release Accessory B			FOLD DOWN TRANSFER ARM SPECIFIC				
CAL	CAL TB133	See pg 9	116	II. C	C1 1			
Z	Multi Cover Use Pricing fo	or Highest Grade	HC	Hinge Cover	Standard			

¹Includes QR bracket connection installed on recliner

⁶Position is locked in place by motors and can only be adjusted by the control wand.







 $^{^2}Standard\ on\ right\ side\ (facing).$ Replace "R" with "L" for left side (facing) option where available.

 $^{^{3} \}mbox{\it Takes}$ place of Push Bar. Paper sold separately.

⁴Contact customer service for solid surface pricing

⁵Over the lap tray function recommended in seated position.