



# DURABLE® POWER SERIES

# ZG Power Healthcare Recliner

## Select Aesthetic Style (Denotes back & arm design)




(Konkert style pictured)

## Functional Style (Aesthetic style code precedes selected model number below, e.g., KCZP17)


### CZP17 ZG Power Healthcare Recliner

COM	Dillon/Carson	A <sub>5</sub>	B <sub>6</sub>	C <sub>7</sub>	D <sub>8</sub>	E <sub>9</sub>	F <sub>10</sub>	G <sub>11</sub>	H <sub>12</sub>	I <sub>13</sub>
9199	9479	9739	9999	10259	10529	10919	11319	11839	12369	13029


## Standard Features




Chaise Continuous Leg Support



Recline Functions Operated by Intuitive Control Wand



4" Central Locking Caster System



Weightless Zero Gravity Position via Wand Preset

## Dimensions

	Overall			Seat			Back H	Arm H	Kick Out Seat D	Recline D	Distance from Wall
	W	D	H	W	D	H					
RCZP17   KCZP17	33	37	46.5	24	21	22	27.5	26.5	40.5	68	15.5
KMCZP17	33	37	46.5	24	21	22	27.5	26.5	40.5	68	15.5

### General Info

COM:	6 yds
Static Wt Cap:	650 lb
Dynamic Wt Cap:	500 lb
Unit Weight:	174 lb

## Available Options

### ARM CAP OPTIONS

See pg 2 for finishes

U	Urethane	169
J1	Solid Surface - Group 1	619
J2	Solid Surface - Group 2	819
F	Fully Upholstered Arms	Standard

### CASTER OPTIONS

Z-Q3	3" - Dual Wheel Central Locking System	No Charge
Z-Q4	4" - Central Locking System	Standard
Z-Q5	5" - Central Locking System	119

### HEAT/MASSAGE OPTIONS (IEC 60601-1 CERTIFIED)

<i>Wand Controls</i>		
HW	Heat	1139
HMW	Heat & Massage	1369
<i>In-Arm Controls</i>		
H	Heat	849
HM	Heat & Massage	1119

### ADDITIONAL FEATURES & OPTIONS

BBU	Battery Back Up	649
PB	Push Bar	Standard
EG	Protective Edge Guard	Standard
QR	Additional Quick Release Accessory Bracket	69
CAL	CAL TB133	See pg 9
Z	Multi Cover	Use Pricing for Highest Grade

### PATIENT WELLNESS & COMFORT OPTIONS

ABP	Adjustable Back Pad	459
	Textile Grades COM-B	569
	Textile Grades C-I	339
AHP	Adjustable Head Rest Pillow	498
DP	Dual IV Poles & Brackets <sup>1</sup>	89
FBH	Foley Bag Holder	279
FT	Foot Tray	249
IV	IV Pole & Bracket <sup>1</sup>	229
LC	LiquiCell	Standard
LCKFR	Orthopedic Locking Footrest <sup>6</sup>	59
MBB	Moisture Barrier (Back)	59
MBS	Moisture Barrier (Seat)	849
MR	Monitor Arm (right side facing) <sup>1,2</sup>	949
MTR	Monitor Arm with Tray (right side facing) <sup>1,2</sup>	229
PR1	Paper Roll Holder <sup>3</sup>	Standard
TE	Trendelenburg	459
USB	USB Power Outlet (IEC 60601-1 Certified)	229
WR	Wrist Restraint Mounts (pair)	

### TRAY TABLE OPTIONS

TSR1	Adjustable (Right) Side Tray <sup>1,2,5</sup>	399
TD3	Dual Adjustable Side Trays (TSR1 & TSL1) <sup>1,5</sup>	798
LFSTR	Large Folding (Right) Side Tray <sup>4</sup>	249
LFSTRCH	Large Folding Cup Holder (Right) Side Tray	319
PLFSTR	Padded Large Folding (Right) Side Tray	379
SFSTR	Small Folding (Right) Side Tray <sup>4</sup>	249

\*Refer to pg 103 for full option details and images\*

<sup>1</sup>Includes QR bracket connection installed on recliner

<sup>2</sup>Standard on right side (facing). Replace "R" with "L" for left side (facing) option where available.

<sup>3</sup>Takes place of Push Bar. Paper sold separately.

<sup>4</sup>Contact customer service for solid surface pricing

<sup>5</sup>Over the lap tray function recommended in seated position.

<sup>6</sup>Position is locked in place by motors and can only be adjusted by the control wand.