



**Select Aesthetic Style** (Denotes back & arm design)



Dual Fold Down Arms

**Select Functional Style** (Aesthetic style code precedes selected model number below, e.g., KC117)

**4062 Mobile Healthcare Recliner - Narrow Arms - Infinite Position Non-Locking Back - Dual Fold Down Arms**

COM	Dillon/Carson	A <sub>5</sub>	B <sub>6</sub>	C <sub>7</sub>	D <sub>8</sub>	E <sub>9</sub>	F <sub>10</sub>	G <sub>11</sub>	H <sub>12</sub>	I <sub>13</sub>
5,659	5,689	5,769	5,979	6,179	6,369	6,599	6,839	7,199	7,619	8,039

**4162 Mobile Healthcare Recliner - Narrow Arms - Infinite Position Locking Back - Dual Fold Down Arms**


COM	Dillon/Carson	A <sub>5</sub>	B <sub>6</sub>	C <sub>7</sub>	D <sub>8</sub>	E <sub>9</sub>	F <sub>10</sub>	G <sub>11</sub>	H <sub>12</sub>	I <sub>13</sub>
6,369	6,419	6,459	6,629	6,839	7,019	7,239	7,479	7,719	7,969	8,429

**C4062 Mobile Healthcare Chaise Recliner - Narrow Arms - Infinite Position Non-Locking Back - Dual Fold Down Arms**


COM	Dillon/Carson	A <sub>5</sub>	B <sub>6</sub>	C <sub>7</sub>	D <sub>8</sub>	E <sub>9</sub>	F <sub>10</sub>	G <sub>11</sub>	H <sub>12</sub>	I <sub>13</sub>
5,859	5,939	6,079	6,319	6,609	6,859	7,189	7,539	7,899	8,069	8,549

**C4162 Mobile Healthcare Chaise Recliner - Narrow Arms - Infinite Position Locking Back - Dual Fold Down Arms**

COM	Dillon/Carson	A <sub>5</sub>	B <sub>6</sub>	C <sub>7</sub>	D <sub>8</sub>	E <sub>9</sub>	F <sub>10</sub>	G <sub>11</sub>	H <sub>12</sub>	I <sub>13</sub>
6,599	6,629	6,769	7,019	7,279	7,529	7,829	8,169	8,569	8,879	9,349



Chaise Models



Locking Back Lever (located on outside arm on applicable models)

**Dimensions**

E4062   E4162   EC4062   EC4162	Overall			Seat			Back H	Arm H	Kick Out Seat D	Recline D	Distance from Wall
	W	D	H	W	D	H					
	32.5	34.5	47.5	24	21	21.5	28.5	27	39	75.5	20

**General Info**

COM:	5.7 yds
w/ Chaise COM:	6.1 yds
Wt Cap:	650 lb
Unit Weight:	155 lb

**Available Options**

**ARM CAP OPTIONS**

See pg 2 for finishes

U	Urethane	Standard
J1	Solid Surface - Group 1	619
J2	Solid Surface - Group 2	819

**CASTER OPTIONS**

Q3	3" - Dual Wheel Central Locking System	1019
Q4	4" - Central Lock & Steer System	1019
Q5	5" - Central Lock & Steer System	1139
S4	4" - (2) Caster Lock & Steer System	569
4S3R1	4" - (3) Swivel/(1) Rigid Locking Casters	Standard
4S4	4" - (4) Swivel Locking Casters	39

**HEAT/MASSAGE OPTIONS (IEC 60601-1 CERTIFIED)**

Wand Controls

HW	Heat	1139
HMW	Heat & Massage	1369

**ADDITIONAL FEATURES & OPTIONS**

PB	Push Bar	Standard
EG	Protective Edge Guard	Standard
QR	Additional Quick Release Accessory Bracket	69
CAL	CAL TB133	See pg 9
Z	Multi Cover	Use Pricing for Highest Grade

**PATIENT WELLNESS & COMFORT OPTIONS**

ABP	Adjustable Back Pad (COM 1.75 yds)	
	Textile Grades COM-B	459
	Textile Grades C-I	569
AHP	Adjustable Head Rest Pillow (COM 1 yd)	339
DP	Dual IV Poles & Brackets <sup>1</sup>	498
FBH	Foley Bag Holder	89
FT	Foot Tray	279
IV	IV Pole & Bracket <sup>1</sup>	249
LC	LiquiCell	229
LCKFR	Orthopedic Locking Footrest	399
MBB	Moisture Barrier (Back)	59
MBS	Moisture Barrier (Seat)	59
MR67	Monitor Arm, on Fold Down Arm <sup>1,2</sup>	899
MTR67	Monitor Arm with Tray, on Fold Down Arm <sup>1,2</sup>	999
PR1	Paper Roll Holder <sup>3</sup>	229
TE	Trendelenburg	339
USB	USB Power Outlet (IEC 60601-1 Certified)	459
WR	Wrist Restraint Mounts (pair)	229

**TRAY TABLE OPTIONS**

TD362	Dual Adjustable Side Trays, on Fold Down Arms <sup>1,5</sup>	898
TSR167	Adjustable (Right) Side Tray, on Fold Down Arm <sup>1,5</sup>	449
TSL167	Adjustable (Left) Side Tray, on Fold Down Arm <sup>1,5</sup>	449

**FOLD DOWN TRANSFER ARM SPECIFIC**

HC	Hinge Cover	Standard
----	-------------	----------

\*Refer to pg 130 for full option details and images\*

<sup>1</sup>Includes QR bracket connection installed on recliner

<sup>2</sup>Standard on right side (facing). Replace "R" with "L" for left side (facing) option where available.

<sup>3</sup>Takes place of Push Bar. Paper sold separately.

<sup>4</sup>Contact customer service for solid surface pricing

<sup>5</sup>Over the lap tray function recommended in seated position.

