


# DURABLE® 650 SERIES

# Mobile Healthcare Recliners | Patient Transfer

Select Aesthetic Style (Denotes back & arm design)



Single Removable Arm  
(Koncert style pictured)

Select Functional Style (Aesthetic style code precedes selected model number below, e.g., KC117)

**1047 Mobile Healthcare Recliner - Infinite Position Non-Locking Back - Single Removable Arm**

COM	Dillon/Carson	A <sub>5</sub>	B <sub>6</sub>	C <sub>7</sub>	D <sub>8</sub>	E <sub>9</sub>	F <sub>10</sub>	G <sub>11</sub>	H <sub>12</sub>	I <sub>13</sub>
5,019	5,029	5,149	5,329	5,559	5,749	5,999	6,259	6,629	6,929	7,269

**1147 Mobile Healthcare Recliner - Infinite Position Locking Back - Single Removable Arm**


COM	Dillon/Carson	A <sub>5</sub>	B <sub>6</sub>	C <sub>7</sub>	D <sub>8</sub>	E <sub>9</sub>	F <sub>10</sub>	G <sub>11</sub>	H <sub>12</sub>	I <sub>13</sub>
5,299	5,339	5,449	5,639	5,859	6,079	6,329	6,589	6,919	7,269	7,609

**C1047 Mobile Healthcare Chaise Recliner - Infinite Position Non-Locking Back - Single Removable Arm**

COM	Dillon/Carson	A <sub>5</sub>	B <sub>6</sub>	C <sub>7</sub>	D <sub>8</sub>	E <sub>9</sub>	F <sub>10</sub>	G <sub>11</sub>	H <sub>12</sub>	I <sub>13</sub>
5,199	5,289	5,469	5,699	5,999	6,259	6,599	6,949	7,479	7,849	8,299

**C1147 Mobile Healthcare Chaise Recliner - Infinite Position Locking Back - Single Removable Arm**

COM	Dillon/Carson	A <sub>5</sub>	B <sub>6</sub>	C <sub>7</sub>	D <sub>8</sub>	E <sub>9</sub>	F <sub>10</sub>	G <sub>11</sub>	H <sub>12</sub>	I <sub>13</sub>
5,519	5,589	5,729	6,019	6,279	6,579	6,929	7,269	7,749	8,069	8,549



Chaise Models

**NOTE:** Transfer arm standard on right side (facing). If specified on left side, handle will be placed on right side and vice versa.

## Dimensions

E1047   E1147   EC1047   EC1147	Overall			Seat			Back H	Arm H	Kick Out Seat D	Recline D	Distance from Wall
	W	D	H	W	D	H					
R1047   R1147   RC1047   RC1147	33.5	34.5	47.5	24	21	21.5	28.5	28	40	75.5	20
K1047   K1147   KC1047   KC1147	33.5	38.5	47	24	21	21.5	28	28	40	76.5	20

### General Info

COM: 5.3 yds  
w/ Chaise COM: 5.7 yds  
Wt Cap: 650 lb  
Unit Weight: 166 lb

## Available Options

### ARM CAP OPTIONS

See pg 2 for finishes

U	Urethane	Standard
J1	Solid Surface - Group 1	619
J2	Solid Surface - Group 2	819

### CASTER OPTIONS

Q3	3" - Dual Wheel Central Locking System	1019
Q4	4" - Central Lock & Steer System	1019
Q5	5" - Central Lock & Steer System	1139
S4	4" - (2) Caster Lock & Steer System	569
4S3R1	4" - (3) Swivel/(1) Rigid Locking Casters	Standard
4S4	4" - (4) Swivel Locking Casters	39

### HEAT/MESSAGE OPTIONS (IEC 60601-1 CERTIFIED)

#### Wand Controls

HW	Heat	1139
HMW	Heat & Massage	1369

#### In-Arm Controls

H	Heat	849
HM	Heat & Massage	1119

### ADDITIONAL FEATURES & OPTIONS

PB	Push Bar	Standard
EG	Protective Edge Guard	Standard
QR	Additional Quick Release Accessory Bracket	69
CAL	CAL TB133	See pg 9
Z	Multi Cover	Use Pricing for Highest Grade

### PATIENT WELLNESS & COMFORT OPTIONS

ABP	Adjustable Back Pad (COM 1.75 yds)	
	Textile Grades COM-B	459
	Textile Grades C-I	569
AHP	Adjustable Head Rest Pillow (COM 1 yd)	339
DP	Dual IV Poles & Brackets <sup>1</sup>	498
FBH	Foley Bag Holder	89
FT	Foot Tray	279
IV	IV Pole & Bracket <sup>1</sup>	249
LC	LiquiCell	229
LCKFR	Orthopedic Locking Footrest	399
MBB	Moisture Barrier (Back)	59
MBS	Moisture Barrier (Seat)	59
MR	Monitor Arm (right side facing) <sup>1,2</sup>	849
MTR	Monitor Arm with Tray (right side facing) <sup>1,2</sup>	949
PR1	Paper Roll Holder <sup>3</sup>	229
TE	Trendelenburg	339
USB	USB Power Outlet (IEC 60601-1 Certified)	459
WR	Wrist Restraint Mounts (pair)	229

### TRAY TABLE OPTIONS

TSR1	Adjustable (Right) Side Tray <sup>1,2,5</sup>	399
TD3	Dual Adjustable Side Trays (TSR1 & TSL1) <sup>1,5</sup>	798

\*Refer to pg 130 for full option details and images\*

<sup>1</sup>Includes QR bracket connection installed on recliner

<sup>2</sup>Standard on right side (facing). Replace "R" with "L" for left side (facing) option where available.

<sup>3</sup>Takes place of Push Bar. Paper sold separately.

<sup>4</sup>Contact customer service for solid surface pricing

<sup>5</sup>Over the lap tray function recommended in seated position.

